


ADJUSTABLE FLANGE NASAL AIRWAYS

ADJUSTABLE FLANGE RED RUBBER AIRWAY

Description:

The Rusch® Adjustable Flange Red Rubber Airway features:

- Soft red rubber construction
- Adjustable flange for secure placement
- Rounded tip allows for gentle insertion
- Available in a full range of sizes, including pediatric
- Reusable latex product
- Bulk, non-sterile



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Size	Length	Qty/Box
125200-000120	Rusch® Adjustable Flange Red Rubber Airway, 12 Fr	95 mm	5
125200-000140	Rusch® Adjustable Flange Red Rubber Airway, 14 Fr	95 mm	5
125200-000160	Rusch® Adjustable Flange Red Rubber Airway, 16 Fr	95 mm	5
125200-000180	Rusch® Adjustable Flange Red Rubber Airway, 18 Fr	95 mm	5
125200-000200	Rusch® Adjustable Flange Red Rubber Airway, 20 Fr	115 mm	5

Catalog #	Size	Length	Qty/Box
125200-000220	Rusch® Adjustable Flange Red Rubber Airway, 22 Fr	125 mm	5
125200-000240	Rusch® Adjustable Flange Red Rubber Airway, 24 Fr	170 mm	5
125200-000260	Rusch® Adjustable Flange Red Rubber Airway, 26 Fr	170 mm	5
125200-000280	Rusch® Adjustable Flange Red Rubber Airway, 28 Fr	170 mm	5
125200-000300	Rusch® Adjustable Flange Red Rubber Airway, 30 Fr	170 mm	5
125200-000320	Rusch® Adjustable Flange Red Rubber Airway, 32 Fr	170 mm	5
125200-000340	Rusch® Adjustable Flange Red Rubber Airway, 34 Fr	170 mm	5
125200-000360	Rusch® Adjustable Flange Red Rubber Airway, 36 Fr	170 mm	5

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.