

ETCO₂

BI-FLO ETCO₂ SAMPLING CANNULA

Description:

The Bi-Flo ETCO₂ Sampling Cannula allows delivery of supplemental oxygen while simultaneously sampling expired CO₂, featuring a design that allows delivery of oxygen to both nasal prongs and sampling of expired gases from both nares.



To order this product, please call:

+1 (919) 433-4999

Catalog #	Size	Tubing Length	Luer Connector	Qty/Case
1843	Bi-Flo ETCO ₂ Sampling Cannula, Adult	7 ft O ₂ / 3 in. CO ₂	Female	25
1844	Bi-Flo ETCO ₂ Sampling Cannula, Adult	7 ft O ₂ / 7 ft CO ₂	Male	25
1845	Bi-Flo ETCO ₂ Sampling Cannula, Adult	7 ft O ₂ / 7 ft CO ₂	Female	25

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.