


GUEDEL ORAL AIRWAYS

CATH-GUIDE® GUEDEL AIRWAY

Description:

The Hudson RCI® Cath-Guide® Guedel Airway features:

- Flexible vinyl with a reinforced bite block
- Guedel-style with three internal channels
- Individually packaged



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Size	Length	ISO Standard Size Design	Qty/Box
1171	0	55 mm	5.5	48 (Individually Packaged)
1170	1	60 mm	6	48 (Individually Packaged)
1169	2	70 mm	7	48 (Individually Packaged)
1168	3	80 mm	8	48 (Individually Packaged)
1164	4	90 mm	9	48 (Individually Packaged)
1167	5	100 mm	10	48 (Individually Packaged)

Catalog #	Size	Length	ISO Standard Size Design	Qty/Box
1166	6	110 mm	11	48 (Individually Packaged)
1165	7	120 mm	12	48 (Individually Packaged)

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.