

ENDOTRACHEAL TUBE ACCESSORIES

Teleflex offers a variety of endotracheal tube accessories to meet specific patient requirements efficiently and economically.



NO IMAGE
AVAILABLE

HUDSON RCI® ENDOTRACHEAL TUBE HOLDER

The Hudson RCI® Endotracheal Tube Holder helps keep the Hudson RCI® Endotracheal Tube in place without applying tape to the patient's face.

[View Details](#)

NO IMAGE
AVAILABLE

RUSCH® DOUBLE SWIVEL CONNECTOR

The Rusch® Double Swivel Connector is designed with a sealing cap to facilitate the insertion of suction catheters or fiber optic bronchoscopes, and is provided non-sterile.

[View Details](#)

NO IMAGE
AVAILABLE

RUSCH® ENDOTEST™ CUFF PRESSURE MANOMETER

The Rusch® EndoTest™ Cuff Pressure Manometer features:


- One-handed operation
- Accurate, easy-to-read gauge
- Simple finger switch for adjustment

[View Details](#)

NO IMAGE
AVAILABLE

RUSCH® SILKOSPRAY™ UNIVERSAL SILICONE SPRAY

The Rusch® Silkospray™ Universal Silicone Spray reduces friction on the surface of devices made of rubber, latex and PVC. Silkospray is CFC-free, and is not intended for use with devices made of silicone.

[View Details](#)

NO IMAGE
AVAILABLE

SHERIDAN® SHER-I-SWIV/FO™ DOUBLE SWIVEL CONNECTOR

The Sheridan® Sher-I-SWIV/FO™ Double Swivel Connector for use with fiber optics features:

- Special orifice allows easy passage of fiber optic scope
- Maintains seal around scope during bronchoscopy procedure
- Swivel feature reduces breathing circuit torque on the tracheal tube

[View Details](#)