


[Home](#) // [Airway Management](#) // [Airways](#) // [ETCO₂](#)

ETCO₂



SOFTECH® PLUS ETCO₂ SAMPLING CANNULA

The Softech® Plus ETCO₂ Sampling Cannula allows delivery of supplemental oxygen while simultaneously sampling expired CO₂, and features softer nasal prongs and lariat for enhanced patient comfort. The split nare design allows separation of the CO₂ readings and the delivery of oxygen, and helps produce a sharp wave form for clinician diagnostics.

[View Details](#)

SOFTECH® BI-FLO ETCO₂ SAMPLING CANNULA

The Softech® Bi-Flo ETCO₂ Sampling Cannula allows delivery of supplemental oxygen while simultaneously sampling expired CO₂, featuring a design that allows delivery of oxygen to both nasal prongs and sampling of expired gases from both nares.

[View Details](#)

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.