

## HEAT & MOISTURE EXCHANGERS (HME)

### GIBECK® HUMID-FLO® HME

#### Description:

The exclusive design of the Gibeck® Humid-Flo® HME allows it to remain in-line during the first 72 hours of mechanical ventilation, even during aerosol treatments.

- Reduces circuit breaks aiding clinician adherence to clinical guidelines regarding infection control practices
- Protects caregivers from exposure to depressurizing circuit spray
- Allows for the delivery of SVN or MDI aerosols
- Reduces the need to routinely interrupt ventilation, therefore maintains positive-end expiratory pressure
- Hygroscopic bacteriostatic microwell paper media provide optimal moisture output with low resistance to flow
- Individually packaged



To order this product, please call:

**+1 (919) 433-4999**

Catalog #	Description	Qty/Case
19912	Gibeck Humid-Flo HME	25

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.