

HEAT & MOISTURE EXCHANGERS (HME)

GIBECK® TRACH-VENT® HME

Description:

The Gibeck® Trach-Vent® HME and accessory products are designed specifically for tracheostomized or intubated, spontaneously breathing patients.

- Recommended for tidal volumes of 50 mL and above
- For use on both pediatric and adult patients
- Oxy-Vent® Adaptor snaps over Trach-Vent® HME to allow delivery of supplemental oxygen
- Individually packaged



To order this product, please call:

+1 (919) 433-4999

Catalog #	Description	Qty/Case
41112	Trach-Vent HME	50
41131	Trach-Vent T HME	50
41312	Trach-Vent+ HME	50
45512	Oxy-Vent with Tubing	10