

GUEDEL ORAL AIRWAYSNO IMAGE
AVAILABLE**HUDSON RCI® COLOR-CODED GUEDEL AIRWAY**

The Hudson RCI® Color-Coded Guedel Airway features:

- Guedel-style enclosed center channel
- Semi-rigid, nontoxic polyethylene with color-coded bite block
- Individually packaged

[View Details](#)NO IMAGE
AVAILABLE**RUSCH® COLOR-CODED GUEDEL AIRWAY**

The Rusch® Color-Coded Guedel Airway features:

- Color-coded for easy identification
- Smooth finish and rounded edges
- Reinforced bite block helps prevent occlusion

[View Details](#)NO IMAGE
AVAILABLE**SOFT CLEAR GUEDEL AIRWAY**

The Rusch® Soft Clear Guedel Airway features:

- Soft, smooth finish for maximum patient comfort
- Flexible design
- Color-coded bite block

[View Details](#)NO IMAGE
AVAILABLE**CATH-GUIDE® GUEDEL AIRWAY**

The Hudson RCI® Cath-Guide® Guedel Airway features:

- Flexible vinyl with a reinforced bite block
- Guedel-style with three internal channels
- Individually packaged

[View Details](#)NO IMAGE
AVAILABLE**TRADITIONAL GUEDEL AIRWAY**

The Rusch® Traditional Guedel Airway features:

- Rigid plastic
- Economical choice
- Size identified on the flange

[View Details](#)