


RESUSCITATION

LIFESAVER® REUSABLE RESUSCITATOR

Description:

The Hudson RCI® Lifesaver® Reusable Resuscitator is autoclavable (except reservoir bag), easy to assemble, and collapsible for compact, convenient storage. This resuscitator's nonbreathing valve offers low resistance to both inspiratory and expiratory flows. Other features include:

- Adaptor (Item Code 1669) positioned between mask and bag, may be used for pressure or CO₂ monitoring
- Durable masks feature silicone construction
- Adult masks include clear rigid shell



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Description	Qty/Case
5340	Adult (mask and reservoir bag not included)	1
5345	Adult (includes mask, reservoir bag and reservoir valve assembly)	1
5341	Child (mask and reservoir bag not included)	1
5346	Child (includes mask, reservoir bag and reservoir valve assembly)	1
5342	Infant (mask and reservoir bag not included)	1
5347	Infant (includes mask, reservoir bag and reservoir valve assembly)	1

Catalog #	Description	Qty/Case
5336	Expiratory flow diverter	1
51012	Intake valve housing	1
51008	Oxygen reservoir, 600 mL	1
5190	Small Infant Mask (Replacement Part)	1
5191	Large Infant Mask (Replacement Part)	1
5192	Small Child Mask (Replacement Part)	1
5193	Large Child Mask (Replacement Part)	1
5194	Small Adult Mask (Replacement Part)	1
5195	Large Adult Mask (Replacement Part)	1

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.