

LMA® AIRWAY ACCESSORIES

Teleflex provides a selection of LMA® Airway Accessories to facilitate your LMA® Airway procedure.



LMA® CLASSIC EXCEL® ACCESSORIES

LMA® Classic Excel® Airway Removable Connector facilitates fiberoptic assisted intubation.

[View Details](#)

LMA® FASTRACH® REUSABLE AIRWAY ACCESSORIES

The LMA® Fastrach® Reusable Endotracheal Tube (ETT) is designed for use with the LMA Fastrach® Reusable Airway. The LMA® Fastrach® Reusable Stabilizer Rod is used during the removal of the LMA® Fastrach® Reusable Airway after intubation to keep the ETT in place.


- One centimeter markings to assist confirmation of ETT placement following LMA® Fastrach® Airway removal

[View Details](#)

LMA® FASTRACH® SINGLE-USE AIRWAY ACCESSORIES


The LMA® Fastrach® Single-Use Endotracheal Tube (ETT) is designed for use with the LMA® Fastrach® Single-Use Airway. The LMA® Fastrach® Single-Use Stabilizer Rod is used during the removal of the LMA® Fastrach® Single-Use Airway after intubation to keep the ETT in place.

- One centimeter markings to assist confirmation of ETT placement following LMA® Fastrach® Airway removal

[View Details](#)

LMA® FLEXIBLE® AIRWAY ACCESSORIES

The LMA® Flexible® Airway Introducer is designed to facilitate LMA® Flexible® Airway insertion.

[View Details](#)

LMA® PROSEAL® ACCESSORIES

The LMA® ProSeal® Airway Introducer is designed to facilitate insertion of the LMA® ProSeal® Airway.

[View Details](#)