


REUSABLE LMA® AIRWAYS

LMA® CLASSIC EXCEL® AIRWAY

Description:

- First-generation laryngeal mask designed to support fiberscope-assisted direct intubation with endotracheal tubes sized up to 7.5 mm
- Removable connector designed for fiberoptic assisted intubation
- Silicone Cuff - Soft, flexible, and conforms to the anatomy^{1,2}
- Epiglottic Elevating Bar (EEB) designed to aid in intubation
- Reusable up to 60 times³



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Mask Size	Patient Size	Max Cuff Volume (Air)*	Max ETT ID (MM)	Fiberoptic Bronchoscope	Qty/Pack
902030	3	30–50 kg	up to 20 mL	7	5	1
902040	4	50–70 kg	up to 30 mL	7	5	1
902050	5	70–100 kg	up to 40 mL	7.5	5.5	1

*These are maximum volumes that should never be exceeded. It is recommended the cuff be inflated to 60 cm H₂O intracuff pressure. **The LMA Fastrach® Reusable Endotracheal Tube (ETT) has a larger outer diameter than the average ETT—compatible sizes for the LMA Fastrach® Reusable ETT are: LMA Classic Excel® Airway size 3 – 6.5 mm LMA Fastrach® Reusable ETT, LMA Classic Excel® Airway size 4 – 6.5 mm LMA Fastrach® Reusable ETT, LMA Classic Excel® Airway size 5 – 7.0 mm LMA Fastrach® Reusable ETT.

References:

1. Eschertzhuber S, Brimacombe J, Hohlrieder M, Keller C. Anaesthesia. 2009;64(1):79-83.
2. Jagannathan N, Sohn LE, Sawardekar A, Gordon J, Langen KE, Anderson K. Anaesthesia. 2012;67(6):632-639.
3. Teleflex. LMA Classic Excel® Airway Instructions for Use. 2015

regions. Contact customer service to confirm availability in your region.