


## REUSABLE LMA® AIRWAYS

### LMA® CLASSIC® AIRWAY

**Description:**

- General purpose reusable airway for routine and emergency anesthetic procedures
- Silicone Cuff - Soft, flexible, and conforms to the anatomy <sup>1,2</sup>
- Standard straight airway tube
- Tactile inflation pilot balloon
- Reusable up to 40 times<sup>3</sup>



United States

To order this product, please call:

**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**

Lat

Am

Catalog #	Mask Size	Patient Size	Max Cuff Volume (Air)*	Qty/Pack
100010	1	up to 5 kg	up to 4 mL	1
100015	1.5	5–10 kg	up to 7 mL	1
100020	2	10–20 kg	up to 10 mL	1
100025	2.5	20–30 kg	up to 14 mL	1
100030	3	30–50 kg	up to 20 mL	1
100040	4	50–70 kg	up to 30 mL	1

Catalog #	Mask Size	Patient Size	Max Cuff Volume (Air)*	Qty/Pack
100050	5	70–100 kg	up to 40 mL	1
100060	6	over 100 kg	up to 50 mL	1

\*These are maximum volumes that should never be exceeded. It is recommended the cuff be inflated to 60 cm H<sub>2</sub>O intracuff pressure.

**References:**

1. Eschertzhuber S, Brimacombe J, Hohlrieder M, Keller C. *Anaesthesia*. 2009;64(1):79-83.
2. Jagannathan N, Sohn LE, Sawardekar A, Gordon J, Langen KE, Anderson K. *Anaesthesia*. 2012;67(6):632-639.
3. Teleflex. LMA Classic<sup>®</sup> Airway Instructions for Use. 2015

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex<sup>®</sup> item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.