


SINGLE-USE LMA® AIRWAYS

LMA® FASTRACH® AIRWAY

Description:

- Designed as a guide for blind intubation of the trachea without moving the head or neck and allows continuous ventilation between intubation attempts¹
- LMA® Fastrach® Single-Use Airway is temporarily paired with the LMA® Fastrach® Reusable ETT



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Description	Patient Size	Max Cuff Volume (Air)*	Qty/Pack
135130	Size 3 LMA® Fastrach® Single-Use Airway	30 - 50 kg	up to 20 mL	10
135140	Size 4 LMA® Fastrach® Single-Use Airway	50 - 70 kg	up to 30 mL	10
135150	Size 5 LMA® Fastrach® Single-Use Airway	70 - 100 kg	up to 40 mL	10

*These are maximum volumes that should never be exceeded. It is recommended the cuff be inflated to 60 cm H₂O intracuff pressure.

References:

1. Teleflex. LMA Fastrach® Airway Instructions for Use. 2015

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.