


LMA® AIRWAY ACCESSORIES

LMA® FASTRACH® REUSABLE AIRWAY ACCESSORIES

Description:

The LMA® Fastrach® Reusable Endotracheal Tube (ETT) is designed for use with the LMA Fastrach® Reusable Airway. The LMA® Fastrach® Reusable Stabilizer Rod is used during the removal of the LMA® Fastrach® Reusable Airway after intubation to keep the ETT in place.

- One centimeter markings to assist confirmation of ETT placement following LMA® Fastrach® Airway removal



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Description	Max ETT ID (MM)	Qty/Pack
131060	LMA® Fastrach® Reusable ETT	6.0 mm	1
131065	LMA® Fastrach® Reusable ETT	6.5 mm	1
131070	LMA® Fastrach® Reusable ETT	7.0 mm	1
131075	LMA® Fastrach® Reusable ETT	7.5 mm	1
131080	LMA® Fastrach® Reusable ETT	8.0 mm	1
130099	LMA® Fastrach® Reusable Stabilizer Rod	All Sizes	1

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.