


**LMA® AIRWAY ACCESSORIES**

**LMA® FASTRACH® SINGLE-USE AIRWAY ACCESSORIES**

**Description:**

The LMA® Fastrach® Single-Use Endotracheal Tube (ETT) is designed for use with the LMA® Fastrach® Single-Use Airway. The LMA® Fastrach® Single-Use Stabilizer Rod is used during the removal of the LMA® Fastrach® Single-Use Airway after intubation to keep the ETT in place.

- One centimeter markings to assist confirmation of ETT placement following LMA® Fastrach® Airway removal



United States

To order this product, please call:

**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**

Lat  
Am

Catalog #	Description	Max ETT ID (MM)	Qty/Pack
136060	LMA® Fastrach® Single-Use ETT and Stabilizer Rod	6.0 mm	10
136065	LMA® Fastrach® Single-Use ETT and Stabilizer Rod	6.5 mm	10
136070	LMA® Fastrach® Single-Use ETT and Stabilizer Rod	7.0 mm	10
136075	LMA® Fastrach® Single-Use ETT and Stabilizer Rod	7.5 mm	10
136080	LMA® Fastrach® Single-Use ETT and Stabilizer Rod	8.0 mm	10

