


REUSABLE LMA® AIRWAYS

LMA® FLEXIBLE® AIRWAY

Description:

- Wire-reinforced airway tube provides flexibility and resistance to compression
- Ideal for ear, nose, and throat, ophthalmic and dental surgeries, or other procedures requiring shared airway access
- Also available as single-use



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Description	Mask Size	Patient Size	Max Cuff Volume (Air)*	Qty/Pack
110020	LMA Flexible® Reusable Airway	2	10–20 kg	up to 10 mL	1
110025	LMA Flexible® Reusable Airway	2.5	20–30 kg	up to 14 mL	1
110030	LMA Flexible® Reusable Airway	3	30–50 kg	up to 20 mL	1
110040	LMA Flexible® Reusable Airway	4	50–70 kg	up to 30 mL	1
110050	LMA Flexible® Reusable Airway	5	70–100 kg	up to 40 mL	1
110060	LMA Flexible® Reusable Airway	6	over 100 kg	up to 50 mL	1

Catalog #	Description	Mask Size	Patient Size	Max Cuff Volume (Air)*	Qty/Pack
111010	LMA Flexible® Introducer Single-Use Pediatric	2 and 2.5			10
111030	LMA Flexible® Introducer Single-Use Adult	3, 4, 5, and 6			10

*These are maximum volumes that should never be exceeded. It is recommended the cuff be inflated to 60 cm H₂O intracuff pressure.

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.