


REUSABLE LMA® AIRWAYS

LMA® PROSEAL® AIRWAY

Description:

- Versatile reusable LMA® Airway
- Double cuff design enables seal pressures of up to 32 cm H₂O¹
- Second-generation device designed to reduce risk of aspiration with an integrated drain tube and dual seal



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Description	Mask Size	Patient Size	Max Cuff Volume (Air)*	Max Size Orogastic Tube	Qty/Pack
150010	LMA ProSeal® Airway	1	up to 5 kg	up to 4 mL	8 French	1
150015	LMA ProSeal® Airway	1.5	5–10 kg	up to 7 mL	10 French	1
150020	LMA ProSeal® Airway	2	10–20 kg	up to 10 mL	10 French	1
150025	LMA ProSeal® Airway	2.5	20–30 kg	up to 14 mL	14 French	1
150030	LMA ProSeal® Airway	3	30–50 kg	up to 20 mL	16 French	1
150040	LMA ProSeal® Airway	4	50–70 kg	up to 30 mL	16 French	1

Catalog #	Description	Mask Size	Patient Size	Max Cuff Volume (Air)*	Max Size Orogastric Tube	Qty/Pack
150050	LMA ProSeal® Airway	5	70–100 kg	up to 40 mL	18 French	1
151010	LMA ProSeal® Introducer	1-2.5	up to 30 kg			1
151030	LMA ProSeal® Introducer	3-5	30–100 kg			1
15499	LMA ProSeal® Cleaning Brush	all sizes				10

*These are maximum volumes that should never be exceeded. It is recommended the cuff be inflated to 60 cm H₂O intracuff pressure. **LMA® ProSeal® Pro-Pack Kit includes 1 LMA® ProSeal® Airway, 1 LMA® ProSeal® Introducer & 1 LMA® ProSeal® Cuff-Deflator.

References:

1. Cook T.M., Gibbison B. BJA 2007; 99: 436-439.

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.