


SINGLE-USE LMA® AIRWAYS

LMA® UNIQUE EVO™ AIRWAY

Description:

- Available in a full range of sizes
- First-generation laryngeal mask designed to support fiberscope-assisted direct intubation with endotracheal tubes sized up to 8.0 mm
- Available with a standard pilot balloon or with integrated Cuff Pilot™ Technology which is designed to lower risk of post-operative sore throat^{1,2,3}
- Silicone Cuff – Soft, flexible, and conforms to the anatomy^{4,5}
- Fixed-curve airway tube facilitates rapid insertion
- Integrated bite block



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

| Catalog # | Description | Mask Size | Patient Size | Max Intracuff Pressure | Max Endotracheal Tube ID | Qty/Pack |
|-----------|-------------------------|-----------|--------------|------------------------|--------------------------|----------|
| 1D2030 | LMA® Unique EVO™ Airway | 3 | 30-50 kg | 60 cm H ₂ O | 7.0 mm | 10 |
| 1D2040 | LMA® Unique EVO™ Airway | 4 | 50-70 kg | 60 cm H ₂ O | 7.5 mm | 10 |
| 1D2050 | LMA® Unique EVO™ Airway | 5 | 70-100 kg | 60 cm H ₂ O | 8.0 mm | 10 |

It is recommended the cuff to be inflated not exceeding 60 cm H₂O intracuff pressure.

References:

1. Chantzara G, Stroumpoulis K, Alexandrou N, Kokkinos L, Iacovidou N, Xanthos T. *Minerva Anesthesiol.* 2014;80(5):547-555.
2. Seet E, Yousaf F, Gupta S, Subramanyam R, Wong DT, Chung F. *Anesthesiology.* 2010;112(3):652-657.‡
3. Vasanth Karthik R, Ranganathan P, Kulkarni AP, Sharma KS. *J Anesth.* 2014;28(5):662-667.
4. Eschertzhuber S, Brimacombe J, Hohlrieder M, Keller C. *Anaesthesia.* 2009;64(1):79-83.
5. Jagannathan N, Sohn LE, Sawardekar A, Gordon J, Langen KE, Anderson K. *Anaesthesia.* 2012;67(6):632-639.

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.