


## SINGLE-USE LMA® AIRWAYS

### LMA® UNIQUE™ AIRWAY

**Description:**

- Suitable for use in routine or emergency anesthetic procedures, and as a rescue airway during unexpected difficult airway and cardiopulmonary resuscitation situations<sup>1</sup>
- Available in convenience plus pack which includes syringe and lubricant gel
- PVC cuff with standard straight airway tube
- Tactile inflation pilot balloon



United States

To order this product, please call:

**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**

Lat

Am

Catalog #	Description	Mask Size	Patient Size	Max Cuff Volume (Air)*	Qty/Pack
ALAA010SU	LMA® Unique® Single-Use Airway	1	up to 5 kg	up to 4 mL	10
ALAA015SU	LMA® Unique® Single-Use Airway	1.5	5-10 kg	up to 7 mL	10
ALAA020SU	LMA® Unique® Single-Use Airway	2	10-20 kg	up to 10 mL	10
ALAA025SU	LMA® Unique® Single-Use Airway	2.5	20-30 kg	up to 14 mL	10
ALAA030SU	LMA® Unique® Single-Use Airway	3	30-50 kg	up to 20 mL	10
ALAA040SU	LMA® Unique® Single-Use Airway	4	50-70 kg	up to 30 mL	10

Catalog #	Description	Mask Size	Patient Size	Max Cuff Volume (Air)*	Qty/Pack
ALAA050SU	LMA® Unique® Single-Use Airway	5	70-100 kg	up to 40 mL	10
125010	LMA® Unique® Single-Use Airway (Canada only code)	1	Up to 5 kg	up to 4 ml	10
125015	LMA® Unique® Single-Use Airway (Canada only code)	1.5	5-10 kg	up to 7 mL	10
125020	LMA® Unique® Single-Use Airway (Canada only code)	2	10-20 kg	up to 10 mL	10
125025	LMA® Unique® Single-Use Airway (Canada only code)	2.5	20-30 kg	up to 14 mL	10
125030	LMA® Unique® Single-Use Airway (Canada only code)	3	30-50 kg	up to 20 mL	10
125040	LMA® Unique® Single-Use Airway (Canada only code)	4	50-70 kg	up to 30 mL	10
125050	LMA® Unique® Single-Use Airway (Canada only code)	5	70-100 kg	up to 40 mL	10
128010	LMA® Unique® Single-Use Airway Plus Pack (Includes Syringe and Lubricant)	1	up to 5 kg	up to 4 mL	10
128015	LMA® Unique® Single-Use Airway Plus Pack (Includes Syringe and Lubricant)	1.5	5-10 kg	up to 7 mL	10
128020	LMA® Unique® Single-Use Airway Plus Pack (Includes Syringe and Lubricant)	2	10-20 kg	up to 10 mL	10
128025	LMA® Unique® Single-Use Airway Plus Pack (Includes Syringe and Lubricant)	2.5	20-30 kg	up to 14 mL	10
128030	LMA® Unique® Single-Use Airway Plus Pack (Includes Syringe and Lubricant)	3	30-50 kg	up to 20 mL	10
128040	LMA® Unique® Single-Use Airway Plus Pack (Includes Syringe and Lubricant)	4	50-70 kg	up to 30 mL	10
128050	LMA® Unique® Single-Use Airway Plus Pack (Includes Syringe and Lubricant)	5	70-100 kg	up to 40 mL	10

Catalog #	Description	Mask Size	Patient Size	Max Cuff Volume (Air)*	Qty/Pack
12100	LMA® Unique® Airway Pediatric Emergency Kit (Includes Syringe and Lubricant)	1 - 2.5	Up to 30kg	Size Dependent	1 Kit
12200	LMA® Unique® Airway Adult Emergency Kit (Includes Syringe and Lubricant)	3, 4, 5	30-100kg	Size Dependent	1 Kit
12299	LMA® Unique® Airway Complete Emergency Kit (Includes Syringe and Lubricant)	1, 1.5, 2, 2.5	All sizes	Size Dependent	1 Kit

\*These are maximum volumes that should never be exceeded. It is recommended the cuff be inflated to 60 cm H<sub>2</sub>O intracuff pressure. \*\* Syringe and lubricant included in the Plus Pack and Emergency Kit.

**References:**

1. Teleflex. LMA® Unique™ Airway Instructions for Use. 2015.

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.