


ATOMIZATION DEVICES

MAD NASAL™ INTRANASAL MUCOSAL ATOMIZATION DEVICE

Description:

- Facilitates absorption across mucosal membranes to your patients' blood stream without an intravenous line
- The soft, conical plug on the tip forms a seal with the nostril, preventing expulsion of fluid
- The spray atomizes drugs into a fine mist of particles 30-100 microns in size
- Atomizes in any position



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Description	Syringe	Vial Adapter	Qty/Box
MAD100	MAD Nasal™ Intranasal Mucosal Atomization Device with 3 mL Syringe	3 mL	No	25
MAD110	MAD Nasal™ Intranasal Mucosal Atomization Device with 1 mL Syringe	1 mL	No	25
MAD130	MAD Nasal™ Intranasal Mucosal Atomization Device with 1 mL Syringe and Vial Adapter	1 mL	Yes	25
MAD140	MAD Nasal™ Intranasal Mucosal Atomization Device with 3 mL Syringe and Vial Adapter	3 mL	Yes	25
MAD300	MAD Nasal™ Intranasal Mucosal Atomization Device without Syringe	No	No	25

regions. Contact customer service to confirm availability in your region.