


ATOMIZATION DEVICES

MADGIC® LARYNGO-TRACHEAL MUCOSAL ATOMIZATION DEVICE

Description:

- Consistent, fine atomized spray for direct coverage of nose, pharynx, larynx and trachea
- Targeted Coverage - narrow, long, malleable stylet easily reaches vocal cords and conforms to patient anatomy to allow targeted deposition of medication
- Adapts to any luer-lock syringe
- Clinician chooses the medication and exact dose and volume (down to 0.1 mL)



United States
To order this product, please call: 866-246-6990
Canada
To order this product, please call: 800-387-9699
Lat Am

Catalog #	Description	Syringe	Qty/Box
MAD600	MADgic® Laryngo-Tracheal Mucosal Atomization Device with 3 mL Syringe	3 mL	25
MAD700	MADgic® Laryngo-Tracheal Mucosal Atomization Device without Syringe	No	25
MAD720	MADgic® Pediatric Laryngo-Tracheal Mucosal Atomization Device without Syringe	No	25

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.