


ATOMIZATION DEVICES

MADOMIZER® BOTTLE ATOMIZER

Description:

- A simple, effective delivery system for topical anesthetics, vasoconstrictors and other nasal or oral medications
- Unique pump design and disposable applicator tip reduces the risk of patient cross-contamination that can occur with compressed air atomizers
- Design allows targeted delivery of exact drug doses to the nasal and oral mucosa
- Consistent, fine atomized spray (approximately 30 – 100 microns)



United States

To order this product, please call:
866-246-6990

Canada

To order this product, please call:
800-387-9699

Lat
Am

Catalog #	Description	Qty/Box
MAD500	MADomizer® Bottle Atomizer	5
MAD510	MADomizer® Disposable Tip	50
MAD515	MADomizer® Pull Tab Cap	20
MAD520	MADomizer® Replacement Bottle Clip	5