

RESUSCITATION

PEEP VALVES

Description:


Teleflex offers a selection of disposable PEEP Valves for attachment to any resuscitation bag flow diverter or filter with a compatible connection.

30 mm Disposable PEEP Valve (Item Code 5385)

- Adjustable from 5 to 20 cm H₂O
- ISO standard 30 mm ID connector
- Attaches to any resuscitation bag flow diverter or filter with a 30 mm connection

22 mm Disposable PEEP Valve (Item Code 5383)

- Adjustable from 5 to 20 cm H₂O
- ISO standard 22 mm ID connector
- Attaches to any resuscitation bag flow diverter with a 22 mm connection



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Description	Qty/Case
5385	Disposable PEEP Valve – 30 mm	10
5383	Disposable PEEP Valve – 22 mm	10