

ETCO₂

PLUS ETCO₂ SAMPLING CANNULA

Description:

The Plus ETCO₂ Sampling Cannula allows delivery of supplemental oxygen while simultaneously sampling expired CO₂, and features softer nasal prongs and lariat for enhanced patient comfort. The split nare design allows separation of the CO₂ readings and the delivery of oxygen, and helps produce a sharp wave form for clinician diagnostics.



To order this product, please call:
+1 (919) 433-4999

Catalog #	Size	Tubing Length	Luer Connector	Qty/Case
2845	Plus ETCO ₂ Sampling Cannula, Adult	7 ft O ₂ / 7 ft CO ₂	Male	50

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.