

REUSABLE LARYNGOSCOPESNO IMAGE
AVAILABLE**RÜSCH® GREENSPEC™ FIBER OPTIC REUSABLE LARYNGOSCOPE HANDLES**

- Strong, durable stainless steel construction
- Bright Xenon lamp
- Textured to enhance grip

[View Details](#)NO IMAGE
AVAILABLE**RÜSCH® STANDARD REUSABLE BLADES**

- Full range of styles and sizes
- Manufactured with surgical grade stainless steel
- Halogen light

[View Details](#)NO IMAGE
AVAILABLE**RÜSCH® DOLPHIN™ REUSABLE STANDARD HANDLES**

- Water resistant
- Handle can stay intact during disinfection
- Strong, durable stainless steel construction
- Textured to enhance grip

[View Details](#)NO IMAGE
AVAILABLE**RÜSCH® EMERALD® REUSABLE FIBER OPTIC BLADES**

- 5.0 mm diameter fiber optic bundle
- Integrated bundle design for easy cleaning
- Matte finish to reduce glare

[View Details](#)NO IMAGE
AVAILABLE**RÜSCH® REUSABLE STANDARD HANDLES**

- Strong, durable stainless steel construction
- Available in cross cut finish or matte finish handles

[View Details](#)NO IMAGE
AVAILABLE**RÜSCH® SNAPLIGHT® REUSABLE FIBER OPTIC BLADES**

- Snap-on light bundle eliminates screws
- Unique magnetized snap-on light bundle can be removed quickly for easy cleaning
- Integrated oxygen port available on infant sizes

[View Details](#)

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.