

## REUSABLE LMA® AIRWAYS

Teleflex offers a comprehensive line of reusable LMA® Airways, providing anesthesia professionals with a full range of laryngeal mask options, so you have the right airway device for every patient.

NO IMAGE  
AVAILABLE

### LMA® CLASSIC EXCEL® AIRWAY

- First-generation laryngeal mask designed to support fiberscope-assisted direct intubation with endotracheal tubes sized up to 7.5 mm
- Removable connector designed for fiberoptic assisted intubation
- Silicone Cuff - Soft, flexible, and conforms to the anatomy<sup>1,2</sup>
- Epiglottic Elevating Bar (EEB) designed to aid in intubation
- Reusable up to 60 times<sup>3</sup>

[View Details](#)

NO IMAGE  
AVAILABLE

### LMA® CLASSIC® AIRWAY

- General purpose reusable airway for routine and emergency anesthetic procedures
- Silicone Cuff - Soft, flexible, and conforms to the anatomy<sup>1,2</sup>
- Standard straight airway tube
- Tactile inflation pilot balloon
- Reusable up to 40 times<sup>3</sup>

[View Details](#)

NO IMAGE  
AVAILABLE

### LMA® FASTRACH® AIRWAY

- Designed as a guide for blind intubation of the trachea without moving the head or neck and allows continuous ventilation between intubation attempts<sup>1</sup>
- Available in single-use or reusable
- Stabilizer rod and endotracheal tubes included

[View Details](#)

NO IMAGE  
AVAILABLE

### LMA® FLEXIBLE® AIRWAY

- Wire-reinforced airway tube provides flexibility and resistance to compression
- Ideal for ear, nose, and throat, ophthalmic and dental surgeries, or other procedures requiring shared airway access
- Also available as single-use

[View Details](#)

NO IMAGE  
AVAILABLE

### LMA® PROSEAL® AIRWAY

- Versatile reusable LMA® Airway
- Double cuff design enables seal pressures of up to 32 cm H<sub>2</sub>O<sup>1</sup>
- Second-generation device designed to reduce risk of aspiration with an integrated drain tube and dual seal

[View Details](#)