


ENDOBRONCHIAL TUBES

ROBERT SHAW - LEFT BRONCHUS ENDOBRONCHIAL TUBE

Description:

The Rusch® Bronchopart™ Robert Shaw - Left Bronchus Endobronchial Tube features:

- Ruscheliit™ PVC technology provides the right degree of stiffness, flexibility and kink resistance
- Specially shaped bronchial cuff helps prevent accidental right upper lobe occlusion
- High-volume/low-pressure cuffs
- Color-coded for easy port identification, balloons and 15 mm connectors
- Double swivel connections with double caps
- Anatomically-correct curvatures facilitate proper placement
- Single-use, sterile



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Size	Qty/Box
116100260	Rusch® Robert Shaw - Left Bronchus Endobronchial Tube, 26 Fr	1
116100280	Rusch® Robert Shaw - Left Bronchus Endobronchial Tube, 28 Fr	1
116100350	Rusch® Robert Shaw - Left Bronchus Endobronchial Tube, 35 Fr	1
116100370	Rusch® Robert Shaw - Left Bronchus Endobronchial Tube, 37 Fr	1
116100390	Rusch® Robert Shaw - Left Bronchus Endobronchial Tube, 39 Fr	1

Catalog #**Size****Qty/Box**

116100410

Rusch® Robert Shaw - Left Bronchus Endobronchial Tube, 41 Fr

1

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.