


## ENDOBRONCHIAL TUBES

### ROBERT SHAW - RIGHT BRONCHUS ENDOBRONCHIAL TUBE

**Description:**

The Rusch® Bronchopart™ Robert Shaw - Right Bronchus Endobronchial Tube features:

- Ruscheliit™ PVC technology provides the right degree of stiffness, flexibility and kink resistance
- Specially shaped bronchial cuff helps prevent accidental right upper lobe occlusion
- High-volume/low-pressure cuffs
- Color-coded for easy port identification, balloons and 15 mm connectors
- Double swivel connections with double caps
- Anatomically-correct curvatures facilitate proper placement
- Single-use, sterile



United States

To order this product, please call:

**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**

Lat

Am

Catalog #	Size	Qty/Box
116200260	Rusch® Robert Shaw - Right Bronchus Endobronchial Tube, 28 Fr	1
116200350	Rusch® Robert Shaw - Right Bronchus Endobronchial Tube, 35 Fr	1
116200370	Rusch® Robert Shaw - Right Bronchus Endobronchial Tube, 37 Fr	1
116200390	Rusch® Robert Shaw - Right Bronchus Endobronchial Tube, 39 Fr	1
116200410	Rusch® Robert Shaw - Right Bronchus Endobronchial Tube, 41 Fr	1

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.