


SPECIALTY INTUBATION

RUSCH® AGT NASAL PREFORMED CUFFED ETT

Description:

The Rusch® AGT Nasal Preformed Cuffed Endotracheal Tube (ETT) features:

- Incorporates a preformed curvature designed to improve surgical access by directing the tube toward the patient’s forehead
- Preformed shape directs tube over patient’s forehead, reducing pressure on nares
- Thermosensitive Ruschelit™ PVC adapts to the anatomical situation at body temperature
- Single-use, sterile



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat
Am

Catalog #	Size	Qty/Box
111781035	Rusch® AGT Nasal Preformed Cuffed ETT, Size 3.5	10
111781040	Rusch® AGT Nasal Preformed Cuffed ETT, Size 4.0	10
111781045	Rusch® AGT Nasal Preformed Cuffed ETT, Size 4.5	10
111781050	Rusch® AGT Nasal Preformed Cuffed ETT, Size 5.0	10
111781055	Rusch® AGT Nasal Preformed Cuffed ETT, Size 5.5	10
111781060	Rusch® AGT Nasal Preformed Cuffed ETT, Size 6.0	10

Catalog #	Size	Qty/Box
111781065	Rusch® AGT Nasal Preformed Cuffed ETT, Size 6.5	10
111781070	Rusch® AGT Nasal Preformed Cuffed ETT, Size 7.0	10
111781075	Rusch® AGT Nasal Preformed Cuffed ETT, Size 7.5	10
111781080	Rusch® AGT Nasal Preformed Cuffed ETT, Size 8.0	10

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.