


## SPECIALTY INTUBATION

### RUSCH® AGT NASAL PREFORMED UNCUFFED ETT

**Description:**

The Rusch® AGT Nasal Preformed Uncuffed Endotracheal Tube (ETT) features:

- Incorporates a preformed curvature designed to improve surgical access by directing the tube toward the patient's forehead
- Preformed shape directs tube over patient's forehead, reducing pressure on nares
- Thermosensitive Ruschelit™ PVC adapts to the anatomical situation at body temperature
- Single-use, sterile



United States

To order this product, please call:

**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**

Lat  
Am

Catalog #	Size	Qty/Box
100181030	Rusch® AGT Nasal Preformed Uncuffed ETT, Size 3.0	10
100181040	Rusch® AGT Nasal Preformed Uncuffed ETT, Size 4.0	10
100181045	Rusch® AGT Nasal Preformed Uncuffed ETT, Size 4.5	10
100181050	Rusch® AGT Nasal Preformed Uncuffed ETT, Size 5.0	10
100181055	Rusch® AGT Nasal Preformed Uncuffed ETT, Size 5.5	10
100181060	Rusch® AGT Nasal Preformed Uncuffed ETT, Size 6.0	10

<b>Catalog #</b>	<b>Size</b>	<b>Qty/Box</b>
100181065	Rusch® AGT Nasal Preformed Uncuffed ETT, Size 6.5	10

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.