


SPECIALTY INTUBATION

RUSCH® AGT ORAL PREFORMED UNCUFFED ETT

Description:

The Rusch® AGT Oral Preformed Uncuffed Endotracheal Tube (ETT) features:

- Incorporates a preformed curvature designed to reduce tube kinking and improve surgical access by directing the tube toward the patient's chin
- Thermosensitive Ruschelit™ PVC adapts to the anatomical situation at body temperature
- Single-use, sterile



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Size	Qty/Box
100180030	Rusch® AGT Oral Preformed Uncuffed ETT, Size 3.0	10
100180035	Rusch® AGT Oral Preformed Uncuffed ETT, Size 3.5	10
100180040	Rusch® AGT Oral Preformed Uncuffed ETT, Size 4.0	10
100180045	Rusch® AGT Oral Preformed Uncuffed ETT, Size 4.5	10
100180050	Rusch® AGT Oral Preformed Uncuffed ETT, Size 5.0	10
100180055	Rusch® AGT Oral Preformed Uncuffed ETT, Size 5.5	10

Catalog #	Size	Qty/Box
100180060	Rusch® AGT Oral Preformed Uncuffed ETT, Size 6.0	10
100180065	Rusch® AGT Oral Preformed Uncuffed ETT, Size 6.5	10
100180070	Rusch® AGT Oral Preformed Uncuffed ETT, Size 7.0	10

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.