


ENDOTRACHEAL TUBE ACCESSORIES

RUSCH® DOUBLE SWIVEL CONNECTOR

Description:

The Rusch® Double Swivel Connector is designed with a sealing cap to facilitate the insertion of suction catheters or fiber optic bronchoscopes, and is provided non-sterile.



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Description	Qty/Box
514900000	Rusch® Double Swivel Connector, 15 mm	25