


REUSABLE LARYNGOSCOPES

RÜSCH® EMERALD® REUSABLE FIBER OPTIC BLADES

Description:

- 5.0 mm diameter fiber optic bundle
- Integrated bundle design for easy cleaning
- Matte finish to reduce glare



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Type	Size	Fuctional Height	Overall Length	Distal Width	Qty/Box
4449000	Miller	00 Small Premature	7 mm	67 mm	10 mm	1
4450050	Miller	0 Premature	7.5 mm	78 mm	10 mm	1
4451100	Miller	1 Infant	7.5 mm	103 mm	10 mm	1
4451150	Miller	1.5 Infant	7.5 mm	124 mm	11 mm	1
4452200	Miller	2 Child	8.5 mm	154 mm	13 mm	1
4453300	Miller	3 Medium Adult	8.5 mm	195 mm	12.5 mm	1

Catalog #	Type	Size	Fuctional Height	Overall Length	Distal Width	Qty/Box
4454400	Miller	4 Large Adult	10.5 mm	205 mm	17 mm	1
4431100	Macintosh	1 Infant	18 mm	94.5 mm	10 mm	1
4432200	Macintosh	2 Child	16 mm	111 mm	13 mm	1
4433300	Macintosh	3 Medium Adult	19 mm	133 mm	14 mm	1
4433500	Macintosh	3.5 Medium Adult	15 mm	145 mm	15 mm	1
4434400	Macintosh	4 Large Adult	20 mm	151 mm	13 mm	1

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.