


SPECIALTY INTUBATION

RUSCH® FLEX REINFORCED MURPHY CUFFED ETT

Description:

The Rusch® Flex Reinforced Murphy Cuffed Endotracheal Tube (ETT) features:

- Ruscheliit™ PVC technology provides the right degree of stiffness, flexibility and kink resistance
- Tubes are fitted with a reinforced stainless steel spiral core for kink resistance
- Highly flexible with a thin, strong wall
- Premium soft high-volume low-pressure cuff
- Single-use, sterile



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat
Am

| Catalog # | Size | Qty/Box |
|-----------|--|---------|
| 103902035 | Rusch® Flex Reinforced Murphy Cuffed ETT, Size 3.5 | 5 |
| 103902040 | Rusch® Flex Reinforced Murphy Cuffed ETT, Size 4.0 | 5 |
| 103902045 | Rusch® Flex Reinforced Murphy Cuffed ETT, Size 4.5 | 5 |
| 103902050 | Rusch® Flex Reinforced Murphy Cuffed ETT, Size 5.0 | 5 |
| 103902055 | Rusch® Flex Reinforced Murphy Cuffed ETT, Size 5.5 | 5 |
| 103902060 | Rusch® Flex Reinforced Murphy Cuffed ETT, Size 6.0 | 5 |

| Catalog # | Size | Qty/Box |
|------------------|---|----------------|
| 103902065 | Rusch® Flex Reinforced Murphy Cuffed ETT, Size 6.5 | 5 |
| 103902070 | Rusch® Flex Reinforced Murphy Cuffed ETT, Size 7.0 | 5 |
| 103902075 | Rusch® Flex Reinforced Murphy Cuffed ETT, Size 7.5 | 5 |
| 103902080 | Rusch® Flex Reinforced Murphy Cuffed ETT, Size 8.0 | 5 |
| 103902085 | Rusch® Flex Reinforced Murphy Cuffed ETT, Size 8.5 | 5 |
| 103902090 | Rusch® Flex Reinforced Murphy Cuffed ETT, Size 9.0 | 5 |
| 103902095 | Rusch® Flex Reinforced Murphy Cuffed ETT, Size 9.5 | 5 |
| 103902100 | Rusch® Flex Reinforced Murphy Cuffed ETT, Size 10.0 | 5 |

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.