


STYLETS

RUSCH® FLEXI-SLIP™ STYLET

Description:

The Rusch® Flexi-Slip™ Stylet features:

- Plastic coated malleable metal core
- Soft, retracted metal tip
- Coated for easy insertion and withdrawal
- Color-coded
- Single-use, sterile



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat
Am

Catalog #	O.D. Size	Length	Qty/Box
502503	10 Fr (for tube sizes 4.0-6.5)	44 cm	20
502505	12 Fr (for tube sizes 5.0-10.0)	45 cm	20
502507	14 Fr (for tube sizes 7.0-10.0)	45 cm	20

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.