


SPECIALTY INTUBATION

RUSCH® LARYNGOFLEX® LARYNGECTOMY TUBE

Description:

The Rusch® Laryngoflex® Reinforced Laryngectomy Tube features:

- Anatomically shaped laryngectomy tube constructed with Ruschelit™ PVC which provides the right degree of stiffness, flexibility and kink resistance
- Tubes are fitted with a reinforced stainless steel spiral core for kink resistance
- Flexible with a thin, strong wall
- Premium thin walled low-pressure cuff
- Single-use, sterile



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Size	Qty/Box
121181070	Rusch® Laryngoflex® Reinforced Laryngectomy Tube, Size 7.0	2
121181080	Rusch® Laryngoflex® Reinforced Laryngectomy Tube, Size 8.0	2
121181090	Rusch® Laryngoflex® Reinforced Laryngectomy Tube, Size 9.0	2
121181100	Rusch® Laryngoflex® Reinforced Laryngectomy Tube, Size 10.0	2