


SPECIALTY INTUBATION

RUSCH® LASERTUBUS™ MAGILL CUFFED ETT

Description:

The Rusch® LaserTubus™ Magill Cuffed Laser Resistant Endotracheal Tube (ETT) features:

- Made of soft white rubber reinforced with corrugated copper foil and an absorbent sponge
- Double cuff design offers maximum protection
- Manufactured with natural rubber latex



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Size	Qty/Box
102004-000040	Rusch® Laser Tube Laser Resistant Magill Cuffed ETT, Size 4.0	2
102004-000050	Rusch® Laser Tube Laser Resistant Magill Cuffed ETT, Size 5.0	2
102004-000060	Rusch® Laser Tube Laser Resistant Magill Cuffed ETT, Size 6.0	2
102004-000070	Rusch® Laser Tube Laser Resistant Magill Cuffed ETT, Size 7.0	2
102004-000080	Rusch® Laser Tube Laser Resistant Magill Cuffed ETT, Size 8.0	2