


## RESUSCITATION

### RUSCH® MANUAL RESUSCITATION BAG

#### Description:

The Rusch® Manual Resuscitation Bag comes fully assembled and is designed to provide quality, durability and dependability during emergency situations. Designed for the pre-hospital market, the Manual Resuscitator offers robust construction and is fully collapsible for easy and convenient storage.

- Clear, flexible cushion face mask allows for easy patient visualization
- Textured bag for a sure, comfortable grip
- Non-rebreathing valve that directs exhaled gas away from clinician
- Swivel patient connector 15 mm ID; 22 mm OD
- Device not made with natural rubber latex
- Integrated PEEP connection, no adaptor necessary, PEEP valve available separately



United States

To order this product, please call:

**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**

Lat  
Am

Catalog #	Description	Reservoir	PEEP Valve	Qty/Case
157100300	Adult Resuscitation Bag	Bag		6
157300100	Pediatric Resuscitation Bag	Bag		6
157400100	Infant Resuscitation Bag	Bag		6
157100100	Adult Resuscitation Bag	Tube		6
157300300	Pediatric Resuscitation Bag	Tube		6

<b>Catalog #</b>	<b>Description</b>	<b>Reservoir</b>	<b>PEEP Valve</b>	<b>Qty/Case</b>
157200190	PEEP valve – 22 mm		•	10

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.