


MRI LARYNGOSCOPE KIT

RÜSCH® MRI LARYNGOSCOPE KIT

Description:

- Independently certified to support MRI conditions with a static magnetic field of 3 Tesla or less and highest gradient field of 720 Gauss/cm or less
- Unique design — black aluminum handle and gold-plated blades allow the set to be easily distinguishable
- Lightweight handle and blades
- Specially designed carbon zinc dry batteries



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Description	Qty/Box
5852300	Each Kit Contains: Macintosh Blades 2, 3, 4; Miller Blades 1, 2, 3; 1 Medium Fiber Optic Handle; 2 Battery Packs	1 Kit
5853300	Extra Battery Pack	1

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.