


SINGLE-USE LARYNGOSCOPES

RÜSCH® POLARIS™ SINGLE-USE FIBER OPTIC LARYNGOSCOPE BLADES

Description:

- Designed to withstand the forces of intubation with the weight and feel of a reusable blade
- Low-profile blade design
- Matte finish blade helps reduce light reflection
- 4.7 mm enclosed fiber optic cable for greater light transmission*

*In adult sizes



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Size	Overall Length	Distal Width	Weight (w/o batteries)	Color	Qty/Box
4150110	Mac 1	94.5 mm	9.9 mm	57 g	White	10
4150120	Mac 2	111.0 mm	13.0 mm	66 g	Blue	10
4150130	Mac 3	133.0 mm	14.0 mm	79 g	Yellow	10
4150135	Mac 3.5	145.0 mm	14.5 mm	84 g	Red	10
4150140	Mac 4	149.0 mm	13.0 mm	87 g	Pink	10
4150150	Mac 5	159.0 mm	14.5 mm	96 g	Dark Yellow	10

Catalog #	Size	Overall Length	Distal Width	Weight (w/o batteries)	Color	Qty/Box
4150002	Miller 00	66.5 mm	12.8 mm	39 g	Beige	10
4150001	Miller 0	78.0 mm	12.8 mm	43 g	Violet	10
4150010	Miller 1	103.0 mm	12.8 mm	48 g	Orange	10
4150015	Miller 1.5	124.0 mm	12.8 mm	57 g	Dark Blue	10
4150020	Miller 2	154.0 mm	13.5 mm	70 g	Grey	10
4150030	Miller 3	195.0 mm	13.5 mm	82 g	Green	10
4150040	Miller 4	205.0 mm	19.0 mm	107 g	Lavender	10

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.