


FIXED FLANGE NASAL AIRWAYS

RUSCH® PVC NASAL AIRWAYS

Description:

Rusch® PVC Nasal Airways feature:

- Soft PVC construction
- Inner striations for smooth suction catheter insertion
- Beveled tip for comfortable, gentle insertion
- Available in a full range of sizes, including pediatric
- Emergency kits available – include foil packages of lubricating jelly
- Individually packaged, sterile



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Size	Length	Qty/Box
123312	Rusch® PVC Nasal Airway, 12 Fr	60 mm	10
123314	Rusch® PVC Nasal Airway, 14 Fr	72 mm	10
123316	Rusch® PVC Nasal Airway, 16 Fr	85 mm	10
123318	Rusch® PVC Nasal Airway, 18 Fr	105 mm	10
123320	Rusch® PVC Nasal Airway, 20 Fr	115 mm	10

Catalog #	Size	Length	Qty/Box
123322	Rusch® PVC Nasal Airway, 22 Fr	125 mm	10
123324	Rusch® PVC Nasal Airway, 24 Fr	130 mm	10
123326	Rusch® PVC Nasal Airway, 26 Fr	140 mm	10
123328	Rusch® PVC Nasal Airway, 28 Fr	155 mm	10
123330	Rusch® PVC Nasal Airway, 30 Fr	165 mm	10
123332	Rusch® PVC Nasal Airway, 32 Fr	170 mm	10
123334	Rusch® PVC Nasal Airway, 34 Fr	175 mm	10
123336	Rusch® PVC Nasal Airway, 36 Fr	180 mm	10
1233	Rusch® PVC Nasal Airway Emergency Pack 12 Fr-18 Fr (4 Airways)	60 mm-105 mm (4 Airways)	1 Set
1234	Rusch® PVC Nasal Airway Emergency Pack 20 Fr-36 Fr (9 Airways)	115 mm-180 mm (9 Airways)	1 Set

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.