


AIRWAY RESCUE DEVICES

RUSCH® QUICKTRACH® CRICOTHYROTOMY KIT

Description:

The Rusch® QuickTrach® Emergency Cricothyrotomy Kit features:

- Rescue device for failed airways
- Preassembled and ready to use
- Conical plastic cannula with fixation flange and 15 mm connector
- Removable stainless steel cricothyrotomy needle
- Removable safety stopper
- Sterile, single-use



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Description	Size	Qty/Box
120900020	Rusch® QuickTrach® Emergency Cricothyrotomy Kit	2.0	1
120900040	Rusch® QuickTrach® Emergency Cricothyrotomy Kit	4.0	1

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.