


**ENDOTRACHEAL TUBE ACCESSORIES**

**RUSCH® SILKOSPRAY™ UNIVERSAL SILICONE SPRAY**

**Description:**

The Rusch® Silkospray™ Universal Silicone Spray reduces friction on the surface of devices made of rubber, latex and PVC. Silkospray is CFC-free, and is not intended for use with devices made of silicone.



United States

To order this product, please call:

**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**

Lat

Am

| Catalog #     | Description   | Qty/Box |
|---------------|---|---------|
| 556000-000000 | Rusch® Silkospray™ Universal Silicone Spray, 500 mL | 10      |

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.