


## REUSABLE LARYNGOSCOPES

### RÜSCH® SNAPLIGHT® REUSABLE FIBER OPTIC BLADES

**Description:**

- Snap-on light bundle eliminates screws
- Unique magnetized snap-on light bundle can be removed quickly for easy cleaning
- Integrated oxygen port available on infant sizes



United States

To order this product, please call:

**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**

Lat

Am

Catalog #	Type	Size	Qty/Box
2200000	Miller	00 Small Premature	1
2215000	Miller	0 Premature	1
2216100	Miller	1 Infant	1
2217200	Miller	2 Child	1
2218300	Miller	3 Medium Adult	1
2219400	Miller	4 Large Adult	1

<b>Catalog #</b>	<b>Type</b>	<b>Size</b>	<b>Qty/Box</b>
2201100	Macintosh	1 Infant	1
2202200	Macintosh	2 Child	1
2203300	Macintosh	3 Medium Adult	1
2204400	Macintosh	4 Large Adult	1
4457700	Miller OxiPort	00 Small Premature	1
4458800	Miller OxiPort	0 Premature	1
2516100	Replacement Bundle	Fits Item Code 002216100	2
2517200	Replacement Bundle	Fits Item Code 002217200	2
2518300	Replacement Bundle	Fits Item Code 002218300	2
2519400	Replacement Bundle	Fits Item Code 002219400	2
2501100	Replacement Bundle	Fits Item Code 002201100	2
2502200	Replacement Bundle	Fits Item Code 002202200	2
2503300	Replacement Bundle	Fits Item Code 002203300	2
2504400	Replacement Bundle	Fits Item Code 002204400	2
2500000	Replacement Bundle	Fits Item Code 004457700	2

<b>Catalog #</b>	<b>Type</b>	<b>Size</b>	<b>Qty/Box</b>
2515000	Replacement Bundle	Fits Item Code 004458800	2

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.