


REUSABLE LARYNGOSCOPES

RÜSCH® STANDARD REUSABLE BLADES

Description:

- Full range of styles and sizes
- Manufactured with surgical grade stainless steel
- Halogen light



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Type	Size	Fuctional Height	Overall Length	Distal Width	Qty/Box
008600000	Miller	00 Small Premature	10 mm	65 mm	10 mm	1
008615050	Miller	0 Premature	10 mm	78 mm	11 mm	1
008616100	Miller	1 Infant	10 mm	100 mm	11 mm	1
008617200	Miller	2 Child	11 mm	159 mm	13 mm	1
008618300	Miller	3 Medium Adult	11 mm	195 mm	13 mm	1
008619400	Miller	4 Large Adult	15 mm	205 mm	15 mm	1

Catalog #	Type	Size	Fuctional Height	Overall Length	Distal Width	Qty/Box
008601100	Macintosh	1 Infant	14 mm	88 mm	10 mm	1
008602200	Macintosh	2 Child	18 mm	106 mm	11 mm	1
008603300	Macintosh	3 Medium Adult	22 mm	132 mm	15 mm	1
008603350	Macintosh	3.5 Medium Adult	24 mm	145 mm	15 mm	1
008604400	Macintosh	4 Large Adult	26 mm	159 mm	16 mm	1
001967000	Improved View Macintosh	0 Premature	14 mm	89 mm	11 mm	1
001963300	Improved View Macintosh	3 Medium Adult	17 mm	132 mm	14 mm	1
001963350	Improved View Macintosh	3.5 Medium Adult	15.5 mm	147 mm	14.8 mm	1
001964400	Improved View Macintosh	4 Large Adult	18.5 mm	160 mm	14.5 mm	1
008607100	Phillips	1 Infant	13 mm	114 mm	10 mm	1
008608200	Phillips	2 Child	18 mm	166 mm	14 mm	1
008611100	Wisconsin	1 Infant	12 mm	113 mm	12 mm	1
008612200	Wisconsin	2 Child	15 mm	130 mm	15 mm	1
008613300	Wisconsin	3 Medium Adult	15 mm	159 mm	16 mm	1
001967100	Improved View Macintosh	1 Infant	13.5 mm	91 mm	9.5 mm	1

Catalog #	Type	Size	Fuctional Height	Overall Length	Distal Width	Qty/Box
008614400	Wisconsin	4 Large Adult	16.5 mm	200 mm	18 mm	1
001965200	Improved View Macintosh	2 Child	15 mm	110 mm	12 mm	1

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.