


STANDARD INTUBATION

SHERIDAN/HVT® MURPHY EYE ETT

Description:

The Sheridan/HVT® Murphy Eye Endotracheal Tube (ETT) is a high-volume, tapered, low pressure cuff ETT, which features:

- HVT Cuff design facilitates tracheal wall seal at low sealing pressures
- Tapered cuff increases visualization of the vocal chords during intubation
- Single-use, sterile



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Size	Qty/Box
5-10310	Sheridan/HVT® Murphy Eye ETT, Size 5.0	10
5-10311	Sheridan/HVT® Murphy Eye ETT, Size 5.5	10
5-10312	Sheridan/HVT® Murphy Eye ETT, Size 6.0	10
5-10313	Sheridan/HVT® Murphy Eye ETT, Size 6.5	10
5-10314	Sheridan/HVT® Murphy Eye ETT, Size 7.0	10
5-10315	Sheridan/HVT® Murphy Eye ETT, Size 7.5	10

Catalog #	Size	Qty/Box
5-10316	Sheridan/HVT® Murphy Eye ETT, Size 8.0	10
5-10317	Sheridan/HVT® Murphy Eye ETT, Size 8.5	10
5-10318	Sheridan/HVT® Murphy Eye ETT, Size 9.0	10
5-10320	Sheridan/HVT® Murphy Eye ETT, Size 10.0	10

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.