


SPECIALTY INTUBATION

SHERIDAN PREFORMED™ NASAL CUFFED ETT

Description:

The Sheridan Preformed™ Nasal Cuffed Endotracheal Tube (ETT) features:

- Preformed shape directs tube over patient’s forehead, reducing pressure on nares
- Bold mark at center of bend indicates distance to distal tip
- Uncuffed tubes have two Murphy eyes for enhanced patient safety
- Single-use, sterile



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Size	Qty/Box
5-22312	Sheridan® Preformed Nasal Cuffed ETT, Size 6.0	10
5-22313	Sheridan® Preformed Nasal Cuffed ETT, Size 6.5	10
5-22314	Sheridan® Preformed Nasal Cuffed ETT, Size 7.0	10
5-22315	Sheridan® Preformed Nasal Cuffed ETT, Size 7.5	10
5-22316	Sheridan® Preformed Nasal Cuffed ETT, Size 8.0	10