


SPECIALTY INTUBATION

SHERIDAN PREFORMED™ NASAL UNCUFFED ETT

Description:

The Sheridan Preformed™ Nasal Uncuffed Endotracheal Tube (ETT) features:

- Preformed shape directs tube over patient's forehead, reducing pressure on nares
- Bold mark at center of bend indicates distance to distal tip
- Uncuffed tubes have two Murphy eyes for enhanced patient safety
- Single-use, sterile



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

| Catalog # | Size | Qty/Box |
|-----------|--|---------|
| 5-22106 | Sheridan® Preformed Nasal Uncuffed ETT, Size 3.0 | 10 |
| 5-22107 | Sheridan® Preformed Nasal Uncuffed ETT, Size 3.5 | 10 |
| 5-22108 | Sheridan® Preformed Nasal Uncuffed ETT, Size 4.0 | 10 |
| 5-22109 | Sheridan® Preformed Nasal Uncuffed ETT, Size 4.5 | 10 |
| 5-22110 | Sheridan® Preformed Nasal Uncuffed ETT, Size 5.0 | 10 |
| 5-22111 | Sheridan® Preformed Nasal Uncuffed ETT, Size 5.5 | 10 |

| Catalog # | Size | Qty/Box |
|------------------|--|----------------|
| 5-22112 | Sheridan® Preformed Nasal Uncuffed ETT, Size 6.0 | 10 |
| 5-22113 | Sheridan® Preformed Nasal Uncuffed ETT, Size 6.5 | 10 |
| 5-22114 | Sheridan® Preformed Nasal Uncuffed ETT, Size 7.0 | 10 |

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.