


## SPECIALTY INTUBATION

### SHERIDAN PREFORMED™ ORAL CUFFED ETT

**Description:**

The Sheridan Preformed™ Oral Cuffed Endotracheal Tube (ETT) features:

- Preformed shape directs tube downward to rest on patient’s chin without kinking
- Bold mark at center of bend indicates distance to distal tip
- Single-use, sterile



United States

To order this product, please call:

**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**

Lat

Am

Catalog #	Size	Qty/Box
5-22208	Sheridan® Preformed Oral Cuffed ETT, Size 4.0	10
5-22209	Sheridan® Preformed Oral Cuffed ETT, Size 4.5	10
5-22210	Sheridan® Preformed Oral Cuffed ETT, Size 5.0	10
5-22211	Sheridan® Preformed Oral Cuffed ETT, Size 5.5	10
5-22212	Sheridan® Preformed Oral Cuffed ETT, Size 6.0	10
5-22213	Sheridan® Preformed Oral Cuffed ETT, Size 6.5	10

<b>Catalog #</b>	<b>Size</b>	<b>Qty/Box</b>
5-22214	Sheridan® Preformed Oral Cuffed ETT, Size 7.0	10
5-22215	Sheridan® Preformed Oral Cuffed ETT, Size 7.5	10
5-22216	Sheridan® Preformed Oral Cuffed ETT, Size 8.0	10
5-22217	Sheridan® Preformed Oral Cuffed ETT, Size 8.5	10
5-22218	Sheridan® Preformed Oral Cuffed ETT, Size 9.0	10

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.