


SPECIALTY INTUBATION

SHERIDAN PREFORMED™ ORAL UNCUFFED ETT

Description:

The Sheridan Preformed™ Oral Uncuffed Endotracheal Tube (ETT) features:

- Preformed shape directs tube downward to rest on patient's chin without kinking
- Bold mark at center of bend indicates distance to distal tip
- Single-use, sterile



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Size	Qty/Box
5-22006	Sheridan® Preformed Oral Uncuffed ETT, Size 3.0	10
5-22007	Sheridan® Preformed Oral Uncuffed ETT, Size 3.5	10
5-22008	Sheridan® Preformed Oral Uncuffed ETT, Size 4.0	10
5-22009	Sheridan® Preformed Oral Uncuffed ETT, Size 4.5	10
5-22010	Sheridan® Preformed Oral Uncuffed ETT, Size 5.0	10
5-22011	Sheridan® Preformed Oral Uncuffed ETT, Size 5.5	10

Catalog #	Size	Qty/Box
5-22012	Sheridan® Preformed Oral Uncuffed ETT, Size 6.0	10
5-22013	Sheridan® Preformed Oral Uncuffed ETT, Size 6.5	10
5-22014	Sheridan® Preformed Oral Uncuffed ETT, Size 7.0	10

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.